

CHECK LIST FOR PREVENTIVE MAINTENANCE

Checklist for Preventive Maintenance

A. External

1. The exterior is clean
2. It is firm on the floor
3. It is properly leveled
4. Its sides are atleast 10 cm away from walls
5. It is away from direct sunlight
6. Room is well ventilated
7. It is opened only when necessary

B. Internal

1. Lid seals properly without gap
2. Lid seal is clean
3. Ice lining tubes/ice packs are in proper position
4. Ice lining tubes/ice packs filled water to proper level
5. Thickness of frost formation is less than 6 mm
6. Vaccines are neatly placed with space for air circulation
7. DPT, TT and DT are kept in the basket in the ILR and not touching the cooling surface
8. Thermometer has been kept among the vaccine
9. Temperature is recorded twice a day

C. Technical

1. Temperature within prescribed limit (if not, set the thermostat)
2. Voltage stabilizer is working properly and equipment are connected through it.
3. Plug of the voltage stabilizer is fitted properly to the power line
4. Connection of equipment to voltage stabilizer are not lose
5. There is no abnormal noise
6. Compressor mounting bolts are tight